

ICP-MS



sapphire

*Collision Cell Multi-Collector
ICP-MS*



www.nu-ins.com

Sapphire: Technologies Collide for Unparalleled Performance

Sapphire: Collision Cell Multi Collector Inductively Coupled Plasma Mass Spectrometer



Introduction

The application of non-traditional stable isotope systems such as Mg, Si, S, K, Ca, Ti, V, Cr, Fe, Ni and Se in geochemistry, cosmochemistry and life sciences has been spurred by the advancement of the multi-collector ICP-MS (MC-ICP-MS) over the last two decades. This versatile technique, although widely considered as the standard method for the high precision isotopic analysis of such systems, is limited by the molecular species generated by the argon ICP ion source which cause significant interferences on the mass spectrum of these elements. These molecular species can interfere directly with the atomic ions of the same nominal mass, leading to inaccurate isotope ratio determination. A typical solution is to use the high-resolution capabilities of MC-ICP-MS to resolve or partially resolve the isotopes of interest from their respective interferences, but at the cost of significantly reduced ion transmission. Furthermore, analysis of certain isotopes such as ^{40}Ca is not possible on the traditional MC-ICP-MS instruments, as the ^{40}Ar interference cannot be resolved.

Sapphire is a collision/reaction cell MC-ICP-MS incorporating both a “high energy” ion path (traditional MC-ICP-MS) and a separate “low energy” ion path (collision/reaction cell). The instrument offers an ideal solution for the precise and accurate isotopic analysis of both the traditional isotope systems, and the non-traditional isotope systems mentioned above.

Sapphire has combined the Plasma 3 mass analyser with the patented (GB 2535826) dual path transfer optics offering the benefits of collision/reaction cell without compromising performance for classical geological applications.

The ‘dual path’ design allows the cell to be used to remove isobaric interferences for specific applications and can be fully bypassed when analyzing classical isotope systems that do not require removal of interferences or for systems where the high resolution features are preferred.

Application Areas



planetary



life science



geochemistry

find out more www.nu-ins.com

Key features:

Double focussing high precision isotopic measurement mass spectrometer

Enhanced Sensitivity Interface as standard

Dual Path transfer optics, incorporating:

- Low energy path with collision / reaction cell for interference removal
- High energy path for classical multi-collector performance

Up to 4 independently controlled reaction gas inlets

Fully dry pumped configuration available for lower hydrocarbon backgrounds

Sixteen Faraday detectors, with optional software switchable resistor amplifiers and up to six ion-counting detectors (with SEM / Daly options) giving maximum user flexibility and dynamic range

High abundance deceleration filters available for improved abundance sensitivity

High resolution and pseudo high resolution capabilities

Variable Zoom Optics allow for instantaneous switching between isotopic systems during analysis

sapphire

Innovators in Mass Spectrometry



Collision Cell Multi-Collector ICP-MS



environmental



planetary



life science



geochemistry



material science



nuclear

Nu Instruments Limited

Corporate Headquarters
Unit 74, Clywedog Road South,
Wrexham Industrial Estate,
Wrexham,
LL13 9XS - UK.

Tel: +44 (0)1978 661304
nu.info@ametek.com

Innovators in Mass Spectrometry. We deliver cutting-edge science and technology solutions, and offer our customers unparalleled support and maintenance service through the comprehensive **AMECARE** program.

Corporate Headquarters

Nu Instruments UK

• nu.sales@ametek.com

Regional Contacts

Nu Instruments Europe

• nu.europe-sales@ametek.com

Nu Instruments Americas

• nu.americas-sales@ametek.com

Nu Instruments Asia

• nu.asia-sales@ametek.com

visit us online:
www.nu-ins.com

Sapphire version BRO15.3
All mentioned trademarks are registered by their respective owners.

